

Product datasheet for **RC207697L3V**

DNAJC8 (NM_014280) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNAJC8 (NM_014280) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNAJC8
Synonyms:	HSPC331; SPF31
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014280
ORF Size:	759 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207697).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014280.2
RefSeq Size:	1765 bp
RefSeq ORF:	762 bp
Locus ID:	22826
UniProt ID:	O75937
Cytogenetics:	1p35.3
Domains:	DnaJ
MW:	29.7 kDa


[View online »](#)

Gene Summary: Suppresses polyglutamine (polyQ) aggregation of ATXN3 in neuronal cells
(PubMed:27133716).[UniProtKB/Swiss-Prot Function]