

Product datasheet for **RC207699L4V**

E2F4 (NM_001950) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	E2F4 (NM_001950) Human Tagged ORF Clone Lentiviral Particle
Symbol:	E2F4
Synonyms:	E2F-4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001950
ORF Size:	1239 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207699).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001950.3
RefSeq Size:	2100 bp
RefSeq ORF:	1242 bp
Locus ID:	1874
UniProt ID:	Q16254
Cytogenetics:	16q22.1
Domains:	E2F_TDP
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

Protein Pathways: Cell cycle, TGF-beta signaling pathway

MW: 44 kDa

Gene Summary: The protein encoded by this gene is a member of the E2F family of transcription factors. The E2F family plays a crucial role in the control of cell cycle and action of tumor suppressor proteins and is also a target of the transforming proteins of small DNA tumor viruses. The E2F proteins contain several evolutionally conserved domains found in most members of the family. These domains include a DNA binding domain, a dimerization domain which determines interaction with the differentiation regulated transcription factor proteins (DP), a transactivation domain enriched in acidic amino acids, and a tumor suppressor protein association domain which is embedded within the transactivation domain. This protein binds to all three of the tumor suppressor proteins pRB, p107 and p130, but with higher affinity to the last two. It plays an important role in the suppression of proliferation-associated genes, and its gene mutation and increased expression may be associated with human cancer. [provided by RefSeq, Jul 2008]