

Product datasheet for **RC207705L2V**

AFAP1L2 (NM_001001936) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AFAP1L2 (NM_001001936) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AFAP1L2
Synonyms:	CTB-1144G6.4; KIAA1914; XB130
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001001936
ORF Size:	2454 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207705).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001001936.1
RefSeq Size:	4019 bp
RefSeq ORF:	2457 bp
Locus ID:	84632
UniProt ID:	Q8N4X5
Cytogenetics:	10q25.3
MW:	91.3 kDa



[View online »](#)

Gene Summary:

May play a role in a signaling cascade by enhancing the kinase activity of SRC. Contributes to SRC-regulated transcription activation.[UniProtKB/Swiss-Prot Function]