

## Product datasheet for **RC207711L3V**

### FAU (NM\_001997) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | FAU (NM_001997) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | FAU  |
| Synonyms:                 | asr1; FAU1; Fub1; Fubi; MNSFbeta; RPS30; S30   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_001997  |
| ORF Size:                 | 399 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC207711).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001997.3</a>  |
| RefSeq Size:              | 570 bp   |
| RefSeq ORF:               | 402 bp   |
| Locus ID:                 | 2197   |
| UniProt ID:               | <a href="#">P35544</a>   |
| Cytogenetics:             | 11q13.1  |
| Domains:                  | UBQ, Ribosomal_S30   |
| Protein Families:         | Druggable Genome   |



[View online »](#)

**Protein Pathways:** Ribosome

**MW:** 14.2 kDa

**Gene Summary:** This gene is the cellular homolog of the fox sequence in the Finkel-Biskis-Reilly murine sarcoma virus (FBR-MuSV). It encodes a fusion protein consisting of the ubiquitin-like protein fubi at the N terminus and ribosomal protein S30 at the C terminus. It has been proposed that the fusion protein is post-translationally processed to generate free fubi and free ribosomal protein S30. Fubi is a member of the ubiquitin family, and ribosomal protein S30 belongs to the S30E family of ribosomal proteins. Whereas the function of fubi is currently unknown, ribosomal protein S30 is a component of the 40S subunit of the cytoplasmic ribosome and displays antimicrobial activity. Pseudogenes derived from this gene are present in the genome. Similar to ribosomal protein S30, ribosomal proteins S27a and L40 are synthesized as fusion proteins with ubiquitin. [provided by RefSeq, Nov 2014]