

Product datasheet for **RC207742A1V**

Human ETV7 (NM_016135) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ETV7 (NM_016135) AAV Particle
Tag:	Myc-DDK
Symbol:	ETV7
Synonyms:	TEL-2; TEL2; TELB
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207742 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAGGAGGGAGAATTGGCTATTTCTCCTATAAGCCCTGTGGCAGCCATGCCTCCCCTAGGCACCCACG
TGCAAGCCAGATGTGAAGCTCAAATTAACCTGCTGGGTGAAGGGGGATCTGCAAGCTGCCAGGAAGACT
CCGCATCCAGCCCGCACTGTGGAGCAGGGAGGACGTGCTGCACTGGCTGCGCTGGGCAGAGCAGGAGTAC
TCTCTGCCATGCACCGCGGAGCACGGGTTGAGATGAACGGACGCGCCCTCTGCATCCTACCAAGGACG
ACTTCGGGCACCGTGCGCCAGCTCAGGTGACGTCCTGTATGAGCTGCTCCAGTACATCAAGACCCAGCG
GCGAGCCCTGGTGTGTGGACCCTTTTTTGGAGGGATCTTCAGGCTGAAGACGCCACCCAGCACTCTCCA
GTCCCCCGGAAGAGGTGACTGGCCCTCTCAGATGGACACCCGAAGGGGCCACCTGCTGCAGCCACCAG
ACCCAGGGCTTACCAGCAACTTCGGCCACCTGGATGACCCTGGCTGGCAAGGTGGACCCTGGCAAGGA
GGAGTCCCTCACTTATGTCACTGTGCAGAGCTCGGCTGCAGGACCCAGGGGTCTGTTCTTCCCCGCG
ATGCCGCAGGCCCCATTGACGGCAGGATCGCTGACTGCCGCTGCTGTGGGATTACGTGTATCAGCTGC
TCCTTGATACCGATATGAGCCCTACATCAAGTGGGAAGACAAGGACGCCAAGATCTTCCGAGTTGTGGA
TCCAAATGGGCTCGCCAGACTCTGGGGAATCACAAGAACCGGTGAACATGACCTACGAGAAGATGTCT
CGTGCCCTGCGCCACTATTATAAGCTTAATATCATTAAAGAAGAACCGGGGCAGAACTCCTGTTTCAGAT
TTCTAAAGACTCCGGGAAAGATGGTCCAGGACAAGCACAGCCACCTGGAGCCGCTGGAGAGCCAGGAGCA
GGACAGAATAGAGTTCAAGGACAAGAGGCCAGAAATCTCTCCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207742 protein sequence Red=Cloning site Green=Tags(s) MQEGELAISPISPVAAMPPLGTHVQARCEAQINLLGEGGICKLPGRRLRIQPALWSREDVLHWRWAEQEY SLPCTAEHGFEMNGRALCILTKDDFRHRAPSSGDVLYELLQYIKTQRRALVCGPFFGGIFRLKTPTQHSP VPPEEVTGPSQMDTRRGHLLQPPDPGLTSNFGHLDDPGLARWTPGKEESNLCHCAELGCRTQGVCSFPA MPQAPIDGRIADCRLLDYVYQLLLDTRYEPYIKWEDKDAKIFRVDPNGLARLWGNHKNRVNMTYEKMS RALRHYKLNIIKKEPGQKLLFRFLKTPGKMVQDKHSHLEPLESQEQDRIEFKDKRPEISP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_016135
ORF Size:	1023 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_016135.2 , NP_057219.1
RefSeq Size:	1761 bp
RefSeq ORF:	1026 bp
Locus ID:	51513
UniProt ID:	Q9Y603
Cytogenetics:	6p21.31
MW:	39 kDa