

Product datasheet for **RC207750L2V**

R3HDM1 (NM_015361) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	R3HDM1 (NM_015361) Human Tagged ORF Clone Lentiviral Particle
Symbol:	R3HDM1
Synonyms:	R3HDM
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_015361
ORF Size:	3297 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207750).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015361.2
RefSeq Size:	4666 bp
RefSeq ORF:	3300 bp
Locus ID:	23518
UniProt ID:	Q15032
Cytogenetics:	2q21.3
Domains:	R3H
MW:	120.5 kDa



[View online »](#)