

Product datasheet for **RC207754A1V**

Human TP53TG5 (NM_014477) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TP53TG5 (NM_014477) AAV Particle
Tag:	Myc-DDK
Symbol:	TP53TG5
Synonyms:	C20orf10; CLG01
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207754 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGTCCATCAGCAAAGAAGAGGCCCAAGAACAGCAGGGTTTCCAAGATGCAAGATGAGAACTGCGGG
ACGAGACAGAGCAGCCTGTGAGCAAAGTAATTGAGCGGAACCGTCTGAGAACGGTGTTAAAAAATTGTC
GCTCTTGAAGCTACTCAAGAGCTCAAACCGCCGGATCCAAGAAGTGCATAAGCTGGCCAAAAGGTGTTGG
CATTCACTGCTCAGTGTTCCAAGATTCTCCGCATCTCCTCTGGGGAAAAACAGTGCCTGCAATAAAACAA
AACAAAATAATGAAGAGTTCAGGAGATCGGGTCTCCGAGAAGGAAGTCAAGTCCAAGAAATTAGAGTC
CACAGGGGACCCCTAAGAAAAAAGAGTACAAGGAGTGGAAGTCCCAGGTGCAGTCAGGGATGAGGAACAAG
GAGAAAACGTCATTGGCGGCAATGCCACGGAAGAAAAAGCATATAGAGCCTGAGGTTCCAAGGACATCGA
GGGATGATAGCTTGAACCCTGGAGTCCAGGGGAGGCAACCACTCACCGAGGGCCCCGAGTCATCTTCAT
TAAGCCCTACCGCAATAGAACTCCCATGGGGCACATGAAGCAGCTGGATGTAGCCGACCAGTGGATCTGG
TTTGAGGGGCTGCCCACAGAACTCCACCTCCCGCGCCAGGGTGATGTGCAGATCCTCCACCCTGCGTT
GGGTCAAGCGCCGCTGCACCCGCTTCTGCTCCGCATCACTTGAAATGCCTATGTGGCATCCATACAAGGT
TGATGTGACCTGGACGAGAGCCAGAGGTGCGAGCAGAGGGTGGAGATCGCGCCATCAGCTGAAGGGGAGG
AATGGGTGGCGAAATTCAGAGTCTACAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207754 protein sequence Red=Cloning site Green=Tags(s) MSPSAKKRPKNSRVSKMQDEKL RDETEQPVSKVIERNRLRTVLK NLSLLKLLKSSNRRIQELHKLAKRCW HSLLSVPKILRISSGENSACNKTQ NNEEFQEIGCSEKELKSKKLESTGDPKKKEYKEWKSQVQSGMRNK EKTSLAAMP RKEKHIEPEVPRTSRDDSLNPGVQGRQPLTEGPRVIFIKPYRNRTPMGHMKQLDVADQWIW FEGLPTRIHL PAPRVMCRSSTLRWVKRRCTRFCSASLEMPMWHYPYKVDVTWTRARGASRGWRSRHQLKGR NGWRNSRVYK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_014477
ORF Size:	870 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014477.2</u>
RefSeq Size:	995 bp
RefSeq ORF:	873 bp
Locus ID:	27296
UniProt ID:	<u>Q9Y2B4</u>
Cytogenetics:	20q13.12
MW:	34 kDa