

Product datasheet for **RC207768L3V**

Interferon regulatory factor 9 (IRF9) (NM_006084) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Interferon regulatory factor 9 (IRF9) (NM_006084) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Interferon regulatory factor 9
Synonyms:	IRF-9; ISGF3; ISGF3G; p48
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006084
ORF Size:	1179 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207768).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006084.4
RefSeq Size:	1699 bp
RefSeq ORF:	1182 bp
Locus ID:	10379
UniProt ID:	Q00978
Cytogenetics:	14q12
Domains:	IRF



[View online »](#)

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Jak-STAT signaling pathway

MW: 43.7 kDa

Gene Summary: This gene encodes a member of the interferon regulatory factor (IRF) family, a group of transcription factors with diverse roles, including virus-mediated activation of interferon, and modulation of cell growth, differentiation, apoptosis, and immune system activity. Members of the IRF family are characterized by a conserved N-terminal DNA-binding domain containing tryptophan (W) repeats. Mutations in this gene result in Immunodeficiency 65. [provided by RefSeq, Jul 2020]