

## Product datasheet for **RC207780L2V**

### Peroxiredoxin 6 (PRDX6) (NM\_004905) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Peroxiredoxin 6 (PRDX6) (NM_004905) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Peroxiredoxin 6
Synonyms:	1-Cys; aiPLA2; AOP2; HEL-S-128m; LPCAT-5; NSGPx; p29; PRX
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004905
ORF Size:	672 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207780).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004905.2</a>
RefSeq Size:	1715 bp
RefSeq ORF:	675 bp
Locus ID:	9588
UniProt ID:	<a href="#">P30041</a>
Cytogenetics:	1q25.1
Domains:	AhpC-TSA
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Metabolic pathways, Methane metabolism, Phenylalanine metabolism

**MW:** 25 kDa

**Gene Summary:** The protein encoded by this gene is a member of the thiol-specific antioxidant protein family. This protein is a bifunctional enzyme with two distinct active sites. It is involved in redox regulation of the cell; it can reduce H<sub>2</sub>O<sub>2</sub> and short chain organic, fatty acid, and phospholipid hydroperoxides. It may play a role in the regulation of phospholipid turnover as well as in protection against oxidative injury. [provided by RefSeq, Jul 2008]