

## Product datasheet for **RC207781L4V**

### OSCP (ITGBL1) (NM\_004791) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	OSCP (ITGBL1) (NM_004791) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OSCP
Synonyms:	OSCP; TIED
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004791
ORF Size:	1482 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207781).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004791.1</a>
RefSeq Size:	2535 bp
RefSeq ORF:	1485 bp
Locus ID:	9358
UniProt ID:	<a href="#">O95965</a>
Cytogenetics:	13q33.1
Domains:	EGF
MW:	53.9 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a beta integrin-related protein that is a member of the EGF-like protein family. The encoded protein contains integrin-like cysteine-rich repeats. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012]