

Product datasheet for **RC207798A1V**

Human Granzyme K (GZMK) (NM_002104) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Granzyme K (GZMK) (NM_002104) AAV Particle
Tag:	Myc-DDK
Symbol:	Granzyme K
Synonyms:	TRYP2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207798 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACTAAGTTTTCTTCCTTTCTCTGTTTTCTCCTAATAGTTGGGGCTTATATGACTCATGTGTGTTTCA
ATATGGAAATTATTGGAGGGAAGAAGTGTACCTCATTCCAGGCCATTTATGGCCTCCATCCAGTATGG
CGGACATCACGTTTGTGGAGGTGTTCTGATTGATCCACAGTGGGTGCTGACAGCAGCCCACTGCCAATAT
CGGTTTACCAAAGGCCAGTCTCCCACTGTGGTTTTAGGCGCACACTCTCTCTCAAAGAATGAGGCCTCCA
AACAAACACTGGAGATCAAAAAATTTATACCATTCTCAAGAGTTACATCAGATCCTCAATCAAATGATAT
CATGCTGGTTAAGCTTCAAACAGCCGCAAACTCAATAAACATGTCAAGATGCTCCACATAAGATCCAAA
ACCTCTCTTAGATCTGGAACCAATGCAAGGTTACTGGCTGGGAGCCACCGATCCAGATTCATTAAGAC
CTTCTGACACCCTGCGAGAAGTCACTGTTACTGTCTAAGTCGAAAACCTTGCAACAGCCAAAGTTACTA
CAACGGCGACCCTTTTATCACCAAAGACATGGTCTGTGCAGGAGATGCCAAAGGCCAGAAGGATTCCTGT
AAGGGTGACTCAGGGGGCCCTTGATCTGTAAGGTGTCTCCACGCTATAGTCTCTGGAGGTCAATGAAT
GTGGTGTGCCACAAAGCCTGGAATCTACACCCTGTTAACCAAGAAATACCAGACTTGGATCAAAAGCAA
CCTTGTCCTCGCTCATACAAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207798 protein sequence Red=Cloning site Green=Tags(s) MTKFSSFSLFFLIVGAYMTHVCFNMEIIGGKEVSPHSRPFMASIQYGGHHVCGGVLIDPQWVLTAAHCQY RFTKGQSPTVVLGAHSLSKNEASKQTLEIKKFIPFSRVTS DPQSN DIMLVKLQTA AKLNKHKMLHIRSK TSLRSGTKCKVTGWGATDPDSL RPSDTLREVTVTVL SRKLCNSQSYNGDPFITKDMVCAGDAKGQKDSC KGDSGGPLICKGVFHAIVSGGHECGVATKPGIYTLLTKKYQTWIKSNLVPPHTN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002104
ORF Size:	792 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002104.2</u>
RefSeq Size:	1074 bp
RefSeq ORF:	795 bp
Locus ID:	3003
UniProt ID:	<u>P49863</u>
Cytogenetics:	5q11.2
MW:	28.9 kDa