

Product datasheet for **RC207819A1V**

Human SPEM1 (NM_199339) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SPEM1 (NM_199339) AAV Particle
Tag:	Myc-DDK
Symbol:	SPEM1
Synonyms:	C17orf83
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207819 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCATGGTTGAGCGGCCGAGGCCCGAGTGGGCCTCGTATCACAACGCAACAGCTGCCAGG
ACCTGGGCAACTCTGCTCTGTTGCTGCTGGGCCTCATCATCTGCATTAACTTAGCATCAATATAGTGAC
CCTGCTCTGGAGCCGATTCCGTGGTGTCTTATACCAAGTGTCCATGATACCATTTGTGAGAAAGAAGCT
CCCAAGTCATCATTACTCAGAAAGCAGACCCAGCCCCCTAAGAAGCAGAGTTCTCCTGCAGTCCATCTTC
GGTGCACCATGGACCTGTGATGATGACTGTGTCCCGCCCCCAGCTCACCGCCATCGCCGTCGAGGCTC
TCCCACACGCTGTGCTCACTGCCAGTAGCTTGGGCTCCTGACACTGATGACGAGAAGCCTCATCAGTAC
CCAGCCATCTGCTCCTACCACTGGGATGTCCCTGAGGACTGGGAAGGCTTCCAACACACTCAGGGGACCT
GGGTTCCCTGGTCTCAGGATGCCCCAGAGTCCCTCCCCAGACCATCCGCTTCCAGCCTACCGTAGAGGA
AAGGCCCTCAAAACAGGCATATGGTCCGAGCTGGGCCTAAGGGCCTATGTGTATCCTGTGAACCCCCCA
CCTCCCAGCCCTGAGGCTCCTAGCCACAAGAACGGTGGGGAGGGGGCGGTGCCAGAGGCAGAGGCGGCTC
AGTACCAGCCTGTCCAGCTCCACCTGGGCCAGCAGTCATCCCTGAATTTCCCGGCACCGCTCCTC
AGGCCGAATAGTGTATGATGCCCGGGACATGAGACGGCGGCTTCGGGAAGTACCCGGGAGGTGGAGGCC
CTGTCCGGCTGTACCCCTAGCCTCTGGATCCAGCACTGCCGAGGAGACAAGCAAGAATTGGGTGTACC
GTTCCCTAACTGGGAGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207819 protein sequence Red=Cloning site Green=Tags(s) MAMVERPRPEWASYHNCNSNSCQDLGNSVLLLLGLIICINISINIVTLLWSRFRGVLVQVFHDTICEKEA PKSSLLRKQTQPPKKQSSPAVHLRCTMDPVMMTVSPPPAHRHRRRGSPTRCAHCPVAVAPDTDDEKPHQY PAICSYHWDVPEDWEGFQHTQGTWVPWSQDAPESPPQTIRFQPTVEERPLKTGIWSELGLRAYVYPVNP PPSPEAPSHKNGGEGAVPEAAQYQVPAPTLGPAVIEF SRHRSSGRIVYDARDMRRRLREL TREVEA LSGCYPLASGSSTAETS KNWVYRSLTGR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_199339
ORF Size:	927 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_199339.2</u>
RefSeq Size:	1009 bp
RefSeq ORF:	930 bp
Locus ID:	374768
UniProt ID:	<u>Q8N4L4</u>
Cytogenetics:	17p13.1
MW:	34.8 kDa