

Product datasheet for **RC207843A1V**

Human C11ORF77 (HARBI1) (NM_173811) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human C11ORF77 (HARBI1) (NM_173811) AAV Particle
Tag:	Myc-DDK
Symbol:	C11ORF77
Synonyms:	C11orf77
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207843 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTATACCAATAACAGTGCTTGACTGTGACCTCTTGCTATATGGCCGTGGTCACCGGACATTGGACC
GTTTTAAGCTGGATGATGTGACTGATGAATACTTGATGTCCATGTATGGGTTCCACGGCAGTTCATTTA
TTACTTGGTGGAGCTCTTGGGGGCGAATCTTCTAGGCCTACTCAGCGATCCAGGGCTATTAGCCAGAG
ACACAGGTCCTTGACGATTGGGTTTTATACCTCAGGTTCTTCCAGACTCGGATGGGAGATGCCATTG
GAATCAGTCAGGCGCTATGAGTCGTTGTGTGGCAATGTCACCTGAAGCACTTGTGAAAGGGCCCTACA
GTTTCATTCGCTTTCCAGCTGATGAAGCCTCCATTCAAGGCTCTGAAGGATGAATTCTATGGGTTGGCAGGG
ATGCCAGGGGTGATGGGGTGGTTGACTGTATCCATGTGGCCATCAAGGCACCAAATGCTGAAGACCTCT
CCTATGTGAACCGAAAAGGCCTGCATTCTTAACTGCCTGATGGTGTGTGACATTAGAGGGACACTAAT
GACCGTGGAGACAACTGGCCCGGCAGCCTACAGGACTGTGCTGTGCTGCAGCAGTCTTCCCTCAGTAGT
CAGTTTGAAGCGGGTATGCACAAAGATAGCTGGCTTCTGGGTGACAGTTCCTTCTTCTCGAACCTGGC
TCATGACCCCACTTCACATTCTGAACTCCAGCAGAATATCGCTATAACATGGCCCACTTCTGCAACTCA
CAGTGTGATTGAGAAGACTTTCCGAACCTCTGCTCCGATTCCGCTGCCTGGATGGATCCAAGGGGGCA
CTGCAGTACTCACCAGAGAAATCCAGCCATATCATCTTGGCCTGTTGTGCTCTCCACAACATCTCCCTGG
AGCATGGGATGGATGTTTGGTCTCTCCAATGACAGGACCCATGGAACAGCCCCGGAAGAGGAGTATGA
GCATGGAGTCCCTGGACTTAGAGGCTGACCGTATTCGTGAGGAGCTAATGCTCACTCATTTTATGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207843 protein sequence Red=Cloning site Green=Tags(s) MAIPITVLDCDLLLYGRGHRTLDRFKLDDVTDEYLMMSYGFPRQFIYYLVELLGANLSRPTQRSRAISPE TQVLAALGFYTSQSFQTRMGDAIGISQASMSRCVANVTEALVERASQFIRFPADIASIQALKDEFYGLAG MPGVMGVVDCIHVAIKAPNAEDLSYVNRKGLHSLNCLMVCDIRGTLMTVETNWPGLQDCAVLQQSSLSS QFEAGMHKDSWLLGDSSFFLRTWLMTPHLPETPAEYRYNMAHSATHSVIEKTFRTLCSRFRCLDGSKGA LQYSPEKSSHIILACCVLHNISLEHGMDVWSSPMTGPMEQPPEEEYEHMESLDLEADRIQELMLTHFS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_173811
ORF Size:	1047 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_173811.2
RefSeq Size:	1530 bp
RefSeq ORF:	1050 bp
Locus ID:	283254
UniProt ID:	Q96MB7
Cytogenetics:	11p11.2
MW:	39.1 kDa