

## Product datasheet for **RC207871A1V**

### Human WRCH1 (RHOU) (NM\_021205) AAV Particle

#### Product data:

Product Type:	AAV Particles
Product Name:	Human WRCH1 (RHOU) (NM_021205) AAV Particle
Tag:	Myc-DDK
Symbol:	WRCH1
Synonyms:	ARHU; CDC42L1; G28K; hG28K; WRCH1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207871 representing NM_021205 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**

ATGCCCCGCAGCAGGGGGACCCCGCGTTCCCGACCGCTGCGAGGCGCCTCCGGTGCCGCCGCGTCGGGAGCGCGGTGGACGCGGGGACGCGGGCCTGGGGAGCCGGGGGCGGGGGCGTGCGGGGGTGCCGAGGGGCGCGGGCGTCAAGTGCGTGTGGTCGGCGACGGCGGGTGGGCAAGACGAGCCTGGTGGTGAGCTACACCACCAACGGTACCCACCGAGTACATCCCTACTGCCTTCGACAACCTTCCGCGGTGGTGTCTGTGGATGGCGGCCGTGAGACTCCAACCTGTGACACTGCCGGACAGGATGAATTTGACAAGCTGAGGCCCTCTGCTACACCAACACAGACATCTTCTGCTGTGCTTCAGTGTGAGCCCTCATCCTTCCAGAACGTCAGTGAGAAATGGGTGCCGGAGATTGATGCCACTGTCCCAAAGCCCCATCATCCTAGTTGGAACGCAGTCGGATCTCAGAGAAGATGTCAAAGTCCTCATTGAGTTGGACAAATGCAAAGAAAAGCCAGTGCCTGAAGAGGGGCTAAGCTGTGCGCGAGGAAATCAAAGCCGCTCCTACATCGAGTGTTCAGCCTTGACTCAAAAAAACCTCAAAGAGGTCTTTGATGCAGCCATCGTCGCTGGCATTCAATACTCGGACACTCAGCAACAGCCAAAGAAGTCTAAAAGCAGGACTCCAGATAAAATGAAAACTCTCCAAGTCTGGTGAAGAAGTACTGCTGTTTCGTA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Protein Sequence:</b>	>RC207871 representing NM_021205 Red=Cloning site Green=Tags(s)  MPPQQGDPAFPDRCAPPVPPRRERGGRRGPGEPGGRGRAGGAEGRGVKCVLVGDGAVGKTSLVVSYT TNGYPTEYIPTAFDNFSAVVSVSDGRPVRLQLCDTAGQDEFDKLRPLCYTNTDIFLLCFVSVSPSSFQNV EKWVPEIRCHCPKAPIILVGTQSDLREDVKVLIELDKCKEKPVEEAAKLCAEIKAASYIECSALTQKN LKEVFDAIVAGIQYSDTQQPKKSKSRTPDKMKNLSKSWWKKYCCFV  TRTRPLEQKLISEEDLAANDILDYKDDDDKV
<b>Species:</b>	Human
<b>Serotype:</b>	AAV-2
<b>ACCN:</b>	NM_021205
<b>ORF Size:</b>	774 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_021205.4</u>
<b>RefSeq Size:</b>	4369 bp
<b>RefSeq ORF:</b>	777 bp
<b>Locus ID:</b>	58480
<b>UniProt ID:</b>	<u>Q7L0Q8</u>
<b>Cytogenetics:</b>	1q42.13
<b>MW:</b>	28 kDa