

Product datasheet for **RC207878L1V**

PCBP1 (NM_006196) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PCBP1 (NM_006196) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PCBP1 |
| Synonyms: | HEL-S-85; hnRNP-E1; hnRNP-X; HNRPE1; HNRPX |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_006196 |
| ORF Size: | 1068 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC207878). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_006196.2 |
| RefSeq Size: | 1772 bp |
| RefSeq ORF: | 1071 bp |
| Locus ID: | 5093 |
| UniProt ID: | Q15365 |
| Cytogenetics: | 2p13.3 |
| Domains: | KH |
| Protein Pathways: | Spliceosome |



[View online »](#)

MW: 37.5 kDa

Gene Summary: This intronless gene is thought to have been generated by retrotransposition of a fully processed PCBP-2 mRNA. This gene and PCBP-2 have paralogues (PCBP3 and PCBP4) which are thought to have arisen as a result of duplication events of entire genes. The protein encoded by this gene appears to be multifunctional. It along with PCBP-2 and hnRNPk corresponds to the major cellular poly(rC)-binding protein. It contains three K-homologous (KH) domains which may be involved in RNA binding. This encoded protein together with PCBP-2 also functions as translational coactivators of poliovirus RNA via a sequence-specific interaction with stem-loop IV of the IRES and promote poliovirus RNA replication by binding to its 5'-terminal cloverleaf structure. It has also been implicated in translational control of the 15-lipoxygenase mRNA, human Papillomavirus type 16 L2 mRNA, and hepatitis A virus RNA. The encoded protein is also suggested to play a part in formation of a sequence-specific alpha-globin mRNP complex which is associated with alpha-globin mRNA stability. [provided by RefSeq, Jul 2008]