

## Product datasheet for **RC207883L2V**

### GPR86 (P2RY13) (NM\_176894) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GPR86 (P2RY13) (NM_176894) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR86
Synonyms:	FKSG77; GPCR1; GPR86; GPR94; P2Y13; SP174
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_176894
ORF Size:	999 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207883).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_176894.1</a> , <a href="#">NP_795713.1</a>
RefSeq Size:	2769 bp
RefSeq ORF:	1065 bp
Locus ID:	53829
UniProt ID:	<a href="#">Q9BPV8</a>
Cytogenetics:	3q25.1
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction



[View online »](#)

**MW:** 38.4 kDa

**Gene Summary:** The product of this gene belongs to the family of G-protein coupled receptors. This family has several receptor subtypes with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor is activated by ADP. [provided by RefSeq, Sep 2008]