

## Product datasheet for **RC207893L3V**

### **GALNT2 (NM\_004481) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	GALNT2 (NM_004481) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GALNT2
Synonyms:	CDG2T; GalNAc-T2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004481
ORF Size:	1713 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207893).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004481.2</a>
RefSeq Size:	4484 bp
RefSeq ORF:	1716 bp
Locus ID:	2590
UniProt ID:	<a href="#">Q10471</a>
Cytogenetics:	1q42.13
Domains:	RICIN, Glycos_transf_2
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)

**Protein Pathways:** Metabolic pathways, O-Glycan biosynthesis

**MW:** 64.7 kDa

**Gene Summary:** This gene encodes a member of the glycosyltransferase 2 protein family. Members of this family initiate mucin-type O-glycosylation of peptides in the Golgi apparatus. The encoded protein may be involved in O-linked glycosylation of the immunoglobulin A1 hinge region. This gene may influence triglyceride levels, and may be involved Type 2 diabetes, as well as several types of cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014]