

Product datasheet for **RC207896L3V**

RBM26 (NM_022118) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RBM26 (NM_022118) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RBM26
Synonyms:	ARRS2; C13orf10; PPP1R132; PRO1777; SE70-2; ZC3H17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022118
ORF Size:	2940 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207896).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022118.3
RefSeq Size:	5143 bp
RefSeq ORF:	2943 bp
Locus ID:	64062
UniProt ID:	Q5T8P6
Cytogenetics:	13q31.1
Domains:	RRM
MW:	110.7 kDa


[View online »](#)