

Product datasheet for **RC207900L3V**

ZNF96 (ZSCAN12) (NM_001039643) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF96 (ZSCAN12) (NM_001039643) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF96
Synonyms:	dj29K1.2; KIAA0426; ZFP96; ZFP96, ZNF305, ZNF29K1, ZSCAN12, KIAA0426, dj29K1.2; zinc finger and SCAN domain containing 12; zinc finger protein 96; zinc finger protein 305; ZNF29K1; ZNF96; ZNF305
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001039643
ORF Size:	1812 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207900).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001039643.1 , NP_001034732.1
RefSeq Size:	3931 bp
RefSeq ORF:	1814 bp
Locus ID:	9753
Cytogenetics:	6p22.1
Protein Families:	Transcription Factors
MW:	70.2 kDa



[View online »](#)