

Product datasheet for **RC207901A1V**

Human PFD6 (PFDN6) (NM_014260) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PFD6 (PFDN6) (NM_014260) AAV Particle
Tag: Myc-DDK
Symbol: PFD6
Synonyms: H2-KE2; HKE2; KE-2; PFD6
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC207901 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCGGAGCTGATCCAGAAGAAGCTACAGGGAGAAGTGGAGAAATATCAACAGCTACAGAAGGACTTAA
 GTAAATCCATGTCGGGGAGGCAGAACTTGAAGCACAATAACAGAAAATAATATCGTGAAAGAGGAACT
 GGCCCTGCTGGATGGTCCAACGTGGTCTTTAACTTCTGGGTCGGTCTAGTCAAACAGGAGCTGGGG
 GAGGCTCGGGCCACAGTAGGGAAGAGGCTGGACTATATCACAGCTGAAATTAAGCGATACGAATCCCAGC
 TTCGGGATCTTGAGCGGCAGTCAGAGCAACAGAGGGAGACCCTTGCTCAGCTGCAGCAGGAGTCCAGCG
 GGCCAGGCAGCAAAGGCAGGGGCTCCTGGCAAGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC207901 protein sequence
 Red=Cloning site Green=Tags(s)

MAELIQKQLQGEVEKYQQLQKDLKSKSMGRQKLEAQLTENNIVKEELALLDGSNVVFKLLGPVLVKQELG
 EARATVGKRLDYITAEIKRYESQLRDLERQSEQRETLAQLQEFQRAQAAKAGAPGKA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_014260



[View online »](#)

ORF Size:	387 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014260.2</u>
RefSeq Size:	602 bp
RefSeq ORF:	390 bp
Locus ID:	10471
UniProt ID:	<u>O15212</u>
Cytogenetics:	6p21.32
MW:	14.6 kDa