

Product datasheet for **RC207904L1V**

EHF (NM_012153) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EHF (NM_012153) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EHF
Synonyms:	ESE3; ESE3B; ESEJ
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_012153
ORF Size:	900 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207904).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012153.3
RefSeq Size:	5466 bp
RefSeq ORF:	903 bp
Locus ID:	26298
UniProt ID:	Q9NZC4
Cytogenetics:	11p13
Domains:	ETS, SAM_PNT
Protein Families:	Transcription Factors



[View online »](#)

MW: 34.9 kDa

Gene Summary: This gene encodes a protein that belongs to an ETS transcription factor subfamily characterized by epithelial-specific expression (ESEs). The encoded protein acts as a transcriptional repressor and may be involved in epithelial differentiation and carcinogenesis. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2011]