

Product datasheet for **RC207940**

MSANTD4 (NM_032424) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MSANTD4 (NM_032424) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MSANTD4
Synonyms:	KIAA1826
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC207940 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAGCAGTTGAAAAGAAAAGGAAAAGCAATTTTGTGTTCAAGAACTCAGACCCCTTTGAAAGAAA
 TTACGAAAAGGAAAGAAGTCATTTTCCAGCAGCTCAATACAACAATTAATGTGATGAAGCGAATGGC
 TTGGGAAGAGATTGCACAGTGTGTGAATGCTGTAGGAGAAGGAGAACAGAGGACAGGGACAGAGGTGAAA
 AGAAGGTACCTTGACTGGCGAGCACTTATGAAGAGAAAGAGGATGAAGGCCAACATTAAGCTGGTTGGTT
 CAGGATTTCCCTTCCCTCCTCTGATTTGGATGACTCTCTCACTGAAGAGATAGATGAAAAGATTGGATT
 CCGAAATGATGCAAATTTTACTGGCAAATGTGGCAGATTTAGGGATGCAGGTGGATCCCTTAAGTGA
 GTCAAGGTGGAAGAGGAAGAAAGGGATCCGCAGAGTCCTGAATTTGAAATTGAGGAGGAGGAAGAAATGT
 TGTCATCCGTCATACCAGATTCCAGGAGAGAAAATGAACCTCCCGATTTCCCCCACATTGATGAGTTTTT
 TACCCTTAAGTCAACACCATCTAGATCTGCATATGATGAGCCTCATTTGCTCGTAAATATTGAGAAACAG
 AAAGTACAGTTGGAAAACGACGACTGGATATCGAGGCCGAAAGGCTGCAGGTAGAAAAGGAACGCCTAC
 AAATCGAGAAAGAGAGGCTGCGGCATTTAGACATGGAACATGAGCGGCTTCAGCTAGAGAAGGAGCGGCT
 GCAGATTGAAAGAGAAAAGTTGAGGTTACAGATAGTCAATTCAGAGAAACCGTCCTTGAAAATGAACCT
 GGTCAAGGAGAAAAATCCATGCTTCAACACAGGACATAGAAACAGAGAAGTTAAACCTTGAGCGAGAAC
 GCTTGCAACTGGAAAAGGATAGGCTGCAGTTTTTGAAGTTTGAATCTGAGAAGCTGCAGATTGAAAAGGA
 ACGCTTACAGGTAGAGAAAGACAGACTTCAATTCAGAAAGAAGGACACTTGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC207940 protein sequence
Red=Cloning site Green=Tags(s)


MKQLKRKRKSNFSVQETQILLKEITKRKEVIFSKQLNTTINVMKRMAWEEIIAQCVNAVGEQEORTGTVEVK
 RRYLDWRALMKRRMKANIKLVGGFPLPSSDLDDSLTEEIDEKIGFRNDANFDWQNVADFRDAGGSLTE
 VKVEEERDPQSPFEIEEEEEMLSSVIPDSRRENELPDFPHIDEFFTLNSTPSRSAYDEPHLLVNIEKQ
 KLELEKRRLDIEAERLQVEKERLQIEKERLRHLEMEHERLQLEKERLQIEREKRLQVNSEKPSLENEL
 GOGKESMLOPODIETEKLRERLRLEKDRLOFLKFSESKLQIEKERLOVEKDRLRIQKEGHLQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6509_a07.zip

Restriction Sites: Sgfl-Mlul

Cloning Scheme:



ACCN: NM_032424

ORF Size: 1035 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_032424.3

RefSeq Size: 4095 bp

RefSeq ORF: 1038 bp


Locus ID: 84437

UniProt ID: Q8NCY6


Cytogenetics: 11q22.3

MW: 41.1 kDa


Product images:



Circular map for RC207940



Western blot validation of overexpression lysate (Cat# [LY410120]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC207940 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified MSANTD4 protein (Cat# [TP307940]). The protein was produced from HEK293T cells transfected with MSANTD4 cDNA clone (Cat# RC207940) using MegaTran 2.0 (Cat# [TT210002]).