

Product datasheet for **RC207944L1V**

DTX3L (NM_138287) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DTX3L (NM_138287) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DTX3L
Synonyms:	BBAP; RNF143
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_138287
ORF Size:	2220 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207944).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138287.2 , NP_612144.1
RefSeq Size:	5768 bp
RefSeq ORF:	2223 bp
Locus ID:	151636
UniProt ID:	Q8TDB6
Cytogenetics:	3q21.1
Domains:	RING
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Notch signaling pathway

MW: 83.6 kDa

Gene Summary: DTX3L functions as an E3 ubiquitin ligase (Takeyama et al., 2003 [PubMed 12670957]).
[supplied by OMIM, Nov 2009]