

Product datasheet for **RC207974**

TMEM219 (NM_194280) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM219 (NM_194280) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TMEM219
Synonyms:	IGFBP-3R
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC207974 representing NM_194280 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCAACTGCCAGGCAGGGCACAACCTGCACCTGTGTCTGGCCACCACCCACCTCTGGTCTGTGCCA
CTTTGATCCTGCTGCTCCTTGGCCTCTCTGGCCTGGGCCTTGGCAGCTTCCTCCTCACCCACAGGACTGG
CCTGCGCAGCCCTGACATCCCCAGGACTGGGTCTTTTTTTGAGATCTTTTGGCCAGCTGACCCTGTGT
CCCAGGAATGGGACAGTCACAGGGAAGTGGCGAGGGTCTCACGTCGTGGGCTTGCTGACCACCTGAACT
TCGGAGACGGTCCAGACAGGAACAAGACCCGGACATTCCAGGCCACAGTCCTGGGAAGTCAGATGGGATT
GAAAGGATCTTCTGCAGGACAACCTGGTCTTATCACAGCCAGGGTGACCACAGAAAGGACTGCAGGAACC
TGCCATATATTTTAGTGCTGTTCCAGGAATCCTACCCTCCAGCCAGCCACCCATATCCTGCTCAGAGGAGG
GGGCTGGAAATGCCACCCTGAGCCCTAGAATGGGTGAGGAATGTGTTAGTGTCTGGAGCCATGAAGGCCCT
TGTGCTGACCAAGCTGCTCACCTCGGAGGAGCTGGCTCTGTGTGGCTCCAGGCTGCTGGTCTTGGGCTCC
TTCCTGCTTCTCTGTGGCCTTCTCTGCTGTGTCACTGCTATGTGCTTCCACCCGCGCCGGGAGTCCC
ACTGGTCTAGAACCCGGCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC207974 representing NM_194280
Red=Cloning site Green=Tags(s)

MGNCQAGHNLHLCLAHHPPLVCATLILLLLGLSGLGLGSFLLTHRTGLRSPDIPQDWVSFLRSFGQLTLC
 PRNGTVTKWGRGSHVVGLLTLLNFGDGPDRNKTRTFQATVLGSQMGLKSSAGQLVLITARVTTERTAGT
 CLYFSAVPGILPSSQPPISCSSEEGAGNATLSPRMGEECVSVWSHEGLVLTKLLTSEELALCGSRLLVLGS
 FLLLFCGLLCCVTAMCFHPRRESHWSRTRL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1481_g05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_194280

ORF Size: 720 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_194280.4</u>
RefSeq Size:	1107 bp
RefSeq ORF:	723 bp
Locus ID:	124446
UniProt ID:	<u>Q86XT9</u>
Cytogenetics:	16p11.2
Protein Families:	Druggable Genome, Transmembrane
MW:	25.5 kDa
Gene Summary:	Cell death receptor specific for IGFBP3, may mediate caspase-8-dependent apoptosis upon ligand binding.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC207974