

## Product datasheet for **RC207981L2V**

### COG6 (NM\_020751) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	COG6 (NM_020751) Human Tagged ORF Clone Lentiviral Particle
Symbol:	COG6
Synonyms:	CDG2L; COD2; SHNS
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_020751
ORF Size:	1971 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207981).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020751.1</a> , <a href="#">NP_065802.1</a>
RefSeq Size:	3627 bp
RefSeq ORF:	1974 bp
Locus ID:	57511
UniProt ID:	<a href="#">Q9Y2V7</a>
Cytogenetics:	13q14.11
MW:	73.3 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a subunit of the conserved oligomeric Golgi complex that is required for maintaining normal structure and activity of the Golgi apparatus. The encoded protein is organized with conserved oligomeric Golgi complex components 5, 7 and 8 into a sub-complex referred to as lobe B. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2009]