

Product datasheet for **RC207984L3V**

COMMD7 (NM_053041) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	COMMD7 (NM_053041) Human Tagged ORF Clone Lentiviral Particle
Symbol:	COMMD7
Synonyms:	C20orf92; dj1085F17.3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053041
ORF Size:	600 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207984).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053041.1 , NP_444269.1
RefSeq Size:	1928 bp
RefSeq ORF:	603 bp
Locus ID:	149951
UniProt ID:	Q86VX2
Cytogenetics:	20q11.21
MW:	22.6 kDa



[View online »](#)

Gene Summary:

May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). Associates with the NF-kappa-B complex and suppresses its transcriptional activity (PubMed:15799966).[UniProtKB/Swiss-Prot Function]