

## Product datasheet for **RC208005L4V**

### ZNF280C (NM\_017666) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF280C (NM_017666) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF280C
Synonyms:	SUHW3; ZNF633; ZPET
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017666
ORF Size:	2211 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208005).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017666.2</a> , <a href="#">NP_060136.1</a>
RefSeq Size:	4695 bp
RefSeq ORF:	2214 bp
Locus ID:	55609
UniProt ID:	<a href="#">Q8ND82</a>
Cytogenetics:	Xq26.1
Domains:	zf-C2H2
MW:	83.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the zinc finger domain-containing protein family. This family member contains multiple Cys2-His2(C2H2)-type zinc finger domains, the most common type of zinc finger domain that self-folds to form a beta-beta-alpha structure that binds a zinc ion. [provided by RefSeq, Aug 2011]