

## Product datasheet for **RC208007L1V**

### NADPH oxidase 4 (NOX4) (NM\_016931) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NADPH oxidase 4 (NOX4) (NM_016931) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NADPH oxidase 4
Synonyms:	KOX; KOX-1; RENOX
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016931
ORF Size:	1734 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208007).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016931.2</a>
RefSeq Size:	4459 bp
RefSeq ORF:	1737 bp
Locus ID:	50507
UniProt ID:	<a href="#">Q9NPH5</a>
Cytogenetics:	11q14.3
Protein Families:	Druggable Genome, Transmembrane
MW:	66.9 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a member of the NOX-family of enzymes that functions as the catalytic subunit the NADPH oxidase complex. The encoded protein is localized to non-phagocytic cells where it acts as an oxygen sensor and catalyzes the reduction of molecular oxygen to various reactive oxygen species (ROS). The ROS generated by this protein have been implicated in numerous biological functions including signal transduction, cell differentiation and tumor cell growth. A pseudogene has been identified on the other arm of chromosome 11. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Jan 2009]