

## Product datasheet for **RC208014L3V**

### GDNF Receptor alpha 2 (GFRA2) (NM\_001495) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GDNF Receptor alpha 2 (GFRA2) (NM_001495) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GDNF Receptor alpha 2
Synonyms:	GDNFRB; NRTNR-ALPHA; NTNRA; RETL2; TRNR2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001495
ORF Size:	1392 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208014).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001495.4</a>
RefSeq Size:	3376 bp
RefSeq ORF:	1395 bp
Locus ID:	2675
UniProt ID:	<a href="#">O00451</a>
Cytogenetics:	8p21.3
Domains:	GDNF
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 51.5 kDa

**Gene Summary:** Glial cell line-derived neurotrophic factor (GDNF) and neurturin (NTN) are two structurally related, potent neurotrophic factors that play key roles in the control of neuron survival and differentiation. The protein encoded by this gene is a member of the GDNF receptor family. It is a glycosylphosphatidylinositol(GPI)-linked cell surface receptor for both GDNF and NTN, and mediates activation of the RET tyrosine kinase receptor. This encoded protein acts preferentially as a receptor for NTN compared to its other family member, GDNF family receptor alpha 1. This gene is a candidate gene for RET-associated diseases. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]