

## Product datasheet for **RC208020L4V**

### ABHD1 (NM\_032604) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ABHD1 (NM_032604) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABHD1
Synonyms:	LABH1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032604
ORF Size:	1215 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208020).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032604.2</a>
RefSeq Size:	1471 bp
RefSeq ORF:	1218 bp
Locus ID:	84696
UniProt ID:	<a href="#">Q96SE0</a>
Cytogenetics:	2p23.3
Protein Families:	Transmembrane
MW:	45.2 kDa



[View online »](#)

**Gene Summary:**

This gene is a member of the AB hydrolase superfamily and encodes a protein with an alpha/beta hydrolase fold. This domain is common to a number of hydrolytic enzymes of widely differing phylogenetic origins and catalytic functions. [provided by RefSeq, Jul 2008]