

Product datasheet for **RC208027A1V**

Human HNRPH3 (HNRNPH3) (NM_021644) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HNRPH3 (HNRNPH3) (NM_021644) AAV Particle
Tag:	Myc-DDK
Symbol:	HNRPH3
Synonyms:	2H9; HNRPH3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208027 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGATTGGGTTATGAAACATAATGGTCCAAATGACGCTAGTGATGGGACAGTACGACTTCGTGGACTAC
CATTTGGTTGCAGCAAAGAGGAAATAGTTCAGTTCTTTCAAGGGTTGAAATCGTGCCAAATGGGATAAC
ATTGACGATGGACTACCAGGGGAGAAGCACAGGGGAGGCCTTCGTGCAGTTTGCTTCAAAGGAGATAGCA
GAAATGCTCTGGGAAACACAAGGAAAGAATAGGACACAGGTATATTGAGATCTTCAGAAGTAGCAGGA
GTGAAATCAAAGGATTTTATGATCCACCAAGAAGATTGCTGGGACAGCGACCGGACCATATGATAGACC
AATAGGAGGAAGAGGGGTTATTATGGAGCTGGGCGTGAAGTTATGGAGGTTTTGATGACTATGGTGCC
TATAATAATTACGGCTATGGAATGATGGCTTTGATGACAGAATGAGAGATGGAAGAGGTATGGGAGGAC
ATGGCTATGGTGGAGCTGGTATGCAAGTTCAGGTTTTTCATGGTGGTCATTTCTGACATATGAGAGGGTT
GCCTTTTCGTGCAACTGAAATGACATTGCTAATTTCTCTCACCCTAAATCCAATACGAGTTCATATT
GATATTGGAGCTGATGGCAGAGCCACAGGAGAAGCAGATGTAGAGTTTGTGACACATGAAGATGCAGTAG
CTGCCATGTCTAAAGATAAAAAAATGCAACATCGATATTTGAAGTCTTCTGAATTTCTACTCTCTGG
AGGCGGCTCTGGCATGGGAGGTTCTGGAATGGGAGGCTACGGAAGAGATGGAATGGATAATCAGGAGGC
TATGGATCAGTTGGAAGAATGGGAATGGGAACAATTACAGTGGAGGATATGGTACTCTCTGATGGTTTGG
GTGGTTATGGCGTGGTGGTGGAGGCAGTGGAGGTTACTATGGGCAAGCGGCATGAGTGGAGGTGGATG
CGGTGGGATGTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208027 protein sequence Red=Cloning site Green=Tags(s) MDWVMKHNGPNDASDGTVRLRGLPFGCSKEEIVQFFQGLEIVPNGITLTMDYQGRSTGEAFVQFASKEIA ENALGKHKERIGHRYIEIFRSSRSEIKGFYDPPRLLGQRPGPYDRPIGGRGGYYGAGRGSYGGFDDYGG YNNYGYGNDGFDDMRDGRGMGGHGYGGAGDASSGFHGGHFVHMRGLPFRATENDIANFFSPLNPIRVHI DIGADGRATGEADVEFVTHEDAVAAMSKDKNNMQHRYIELFLNSTPGGGSGMGGSGMGGYGRDGMNQGG YGSVGRMGMGNNYSGGYGTPDGLGGYGRGGGSGGYGQGGMSGGWRGMY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_021644
ORF Size:	993 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_021644.2 , NP_067676.2
RefSeq Size:	2382 bp
RefSeq ORF:	996 bp
Locus ID:	3189
UniProt ID:	P31942
Cytogenetics:	10q21.3
MW:	35.2 kDa