

Product datasheet for **RC208030L3V**

BLAP75 (RMI1) (NM_024945) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BLAP75 (RMI1) (NM_024945) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BLAP75
Synonyms:	BLAP75; C9orf76; FAAP75
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024945
ORF Size:	1875 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208030).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024945.2
RefSeq Size:	3508 bp
RefSeq ORF:	1878 bp
Locus ID:	80010
UniProt ID:	Q9H9A7
Cytogenetics:	9q21.32
MW:	70.6 kDa


[View online »](#)

Gene Summary:

RMI1 is a component of protein complexes that limit DNA crossover formation via the dissolution of double Holliday junctions (Raynard et al., 2006 [PubMed 16595695]).[supplied by OMIM, Mar 2008]