

Product datasheet for **RC208037A1V**

Human RNF24 (NM_007219) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RNF24 (NM_007219) AAV Particle
Tag: Myc-DDK
Symbol: RNF24
Synonyms: G1L
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208037 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCTCGGATTTCCACATTACAACCTCAGGATGCCTAATATTGGATTCCAGAATCTGCCTCTCAACA
TATATATTGTGGTTTTTGGTACTGCTATATTTGTCTTCATCCTTAGTTTACTCTTCTGTTGCTACTTGAT
TAGGCTAAGACATCAAGCACACAAGAATTTATGCCTACAAACAGGTTATATTTAAAGAGAAAGTAAAA
GAATTGAATTTACATGAGCTCTGTGCAGTGTGCCTAGAAGACTTCAAGCCTCGAGATGAGTTGGGGATTT
GCCCATGTAAGCACGCCCTCCACAGAAAGTGCCTTATTAAGTGGCTGGAGTTCTGTAAGTGTGTCCCCT
GTGCAACATGCCAGTTCTACAGCTGGCCAGTTGCACAGTAAGCAGGACCGTGGACCCCTCAGGGGCC
CTTCCTGGGGCAGAGAACATTGTA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208037 protein sequence
Red=Cloning site Green=Tags(s)

MSSDFPHYNFRMPNIGFQNLPLNIYIVVFGTAIFVFILSLLFCCYLIRLRHQAHKEFYAYKQVILKEKVK
ELNLHELCAVCLEDFKPRDELGICPKHAFHRKCLIKWLEVRKVCPLCNMPVLQLAQLHSKQDRGPPQGP
LPGAENIV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_007219
ORF Size:	444 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_007219.2</u>
RefSeq Size:	3374 bp
RefSeq ORF:	447 bp
Locus ID:	11237
UniProt ID:	<u>Q9Y225</u>
Cytogenetics:	20p13
MW:	17.2 kDa