

Product datasheet for **RC208044**

IFT52 (NM_016004) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IFT52 (NM_016004) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	IFT52
Synonyms:	C20orf9; CGI-53; NGD2; NGD5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC208044 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGAGAAAGAGCTGCGGAGCACCATTCTTTTCAATGCCTACAAAAAGGAGATATTTACCACCAACAATG
 GCTACAAATCCATGCAGAAAAAACTTCGGAGTAATTGGAAGATTAGAGCTTAAAGATGAAATCACATC
 TGAGAAGTTAAATGGAGTAAACTGTGGATTACAGCTGGGCCAAGGGAATAATTTACTGCAGCTGAGTTT
 GAAATCCTGAAGAAATATCTTGACACTGGTGGAGATGTCTTTGTGATGCTAGGAGAAGGTGGAGAATCCA
 GATTTGACACCAATATTAACCTTTTACTAGAAGAATATGGAATCATGGTTAATAATGATGCTGTGGTTAG
 AAATGTATATCACAATATTTCCATCCTAAAGAAGCTCTAGTTTCCAGTGGAGTCTTGAACAGGGAAATT
 AGCCGAGCTGCAGAAAGGCTGTGCTGGGATCATTGATGAGGAAAGCAGTGGAACAATGCCAGGCTC
 TCACCTTTGTGTATCCTTTTGGTCCACATTGAGTGTGATGAAACCAGCAGTGGCGGTTCTGTCTACAGG
 TTCTGTCTGCTTCCACTTAACAGACCCATTTTGGCTTTCTATCACTCAAAGAACCAAGGTGGGAAGCTG
 GCAGTGTCTGGTTCATGTACATGTTCACTGATCAATATTTGGACAAAGAAGAAAACAGCAAAATCATGG
 ATGTTGTTTTCCAGTGGCTCACGACAGGAGACATCCACCTAAACCAGATTGATGCTGAGGACCCAGAGAT
 TTCTGACTACATGATGCTGCCCTACACAGCCACCTATCAAAGCGGAATCGAGAGTGTCTCCAGGAGAGT
 GATGAGATCCCAAGGGACTTTACCACCCTCTTCGACCTGTCCATCTTCCAGCTGGATACCACCTCCTTCC
 ACAGCGTCATCGAGGCTCACGAGCAGCTAAATGTGAAACATGAACCACTCCAGCTCATCCAGCCTCAGTT
 TGAGACGCCGCTGCCAACCTTCAGCCTGCGGTTTTCTCCAGTTTCCGGGAGTTACCACCTCCTCCT
 CTGGAGCTATTTGATTTAGATGAAACGTTCTCCTCTGAGAAGGCACGGCTGGCTCAGATTACCAATAAGT
 GTACTGAAGAAGACCTGGAATTTTATGTCAGGAAGTGTGGTATATTCTTGGAGTAACCACTAACTACC
 AAAGGACCAACAGGATGCCAAACATATCCTTGAGCACGTCTTCTCCAAGTGGTGGAGTTCAAGAAATTG
 AACCAGGAACATGACATCGATACAAGTGAACAGCATTCCAGAACAATTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC208044 protein sequence
 Red=Cloning site Green=Tags(s)

MEKELRSTILFNAYKKEIFTTNNGYKSMQKKLRSNWKIQSLKDEITSEKLNGVKLWITAGPREKFTAAEF
 EILKKYLDTGDDVFVMLGEGGESRFDTNINFLLEEYIMVNDAVVRNVYHKYFHPKEALVSSGLNREI
 SRAAGKAVPGIIDEESSGNNALTFVYPFGATLSVMKPAVAVLSTGSVCFLNRPILAFYHSKNQGGKL
 AVLGSCHMFSDQYLDKEENSKIMDVVFQWLTTGDIHLNQIDAEDPEISDYMMPLPYATLSKRNRCLQES
 DEIPRDFTTLFDLSIFQLDTSFHSVIEAHEQLNVKHEPLQLIQPFETPLPTLQPAVFPPSFRELPPPP
 LELFDLDETFSSSEKARLAQITNKCTEEDLEFYVRKCGDILGVTSKLPKDQDQDAKHILEHVFFQVVEFKKL
 NQEHIDTSETAFQNNF


TRTRPLE**QKLISEEDLAANDILDYKDDDDKV**

Chromatograms:

https://cdn.origene.com/chromatograms/mk6088_h02.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_016004

ORF Size: 1311 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_016004.2](#), [NP_057088.2](#)

RefSeq Size: 1675 bp

RefSeq ORF: 1314 bp

Locus ID: 51098


UniProt ID: [Q9Y366](#)

Cytogenetics: 20q13.12


MW: 49.7 kDa

Gene Summary: This gene encodes a conserved proline-rich protein that is a component of the intraflagellar transport-B (IFT-B) core complex. The encoded protein is essential for the integrity of the IFT-B core complex, and for biosynthesis and maintenance of cilia. Mutations in this gene are associated with ciliopathy that affects the skeleton. [provided by RefSeq, Oct 2016]


Product images:



Circular map for RC208044



Western blot validation of overexpression lysate (Cat# [LY402484]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC208044 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified IFT52 protein (Cat# [TP308044]). The protein was produced from HEK293T cells transfected with IFT52 cDNA clone (Cat# RC208044) using MegaTran 2.0 (Cat# [TT210002]).