

Product datasheet for **RC208044L3V**

IFT52 (NM_016004) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IFT52 (NM_016004) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IFT52
Synonyms:	C20orf9; CGI-53; NGD2; NGD5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016004
ORF Size:	1311 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208044).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016004.2 , NP_057088.2
RefSeq Size:	1675 bp
RefSeq ORF:	1314 bp
Locus ID:	51098
UniProt ID:	Q9Y366
Cytogenetics:	20q13.12
MW:	49.7 kDa


[View online »](#)

Gene Summary:

This gene encodes a conserved proline-rich protein that is a component of the intraflagellar transport-B (IFT-B) core complex. The encoded protein is essential for the integrity of the IFT-B core complex, and for biosynthesis and maintenance of cilia. Mutations in this gene are associated with ciliopathy that affects the skeleton. [provided by RefSeq, Oct 2016]