

## Product datasheet for **RC208048L4V**

### HOXA11 (NM\_005523) Human Tagged ORF Clone Lentiviral Particle


#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HOXA11 (NM_005523) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HOXA11
Synonyms:	HOX1; HOX11; RUSAT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005523
ORF Size:	939 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208048).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005523.4</a>
RefSeq Size:	2295 bp
RefSeq ORF:	942 bp
Locus ID:	3207
UniProt ID:	<a href="#">P31270</a>
Cytogenetics:	7p15.2
Protein Families:	Transcription Factors
MW:	34.3 kDa


[View online »](#)

**Gene Summary:**

In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. This gene is involved in the regulation of uterine development and is required for female fertility. Mutations in this gene can cause radio-ulnar synostosis with amegakaryocytic thrombocytopenia. [provided by RefSeq, Jul 2008]

**Product images:**


[RC208048L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC208048L4V particle to overexpress human HOXA11-mGFP fusion protein.