

Product datasheet for **RC208061L3V**

CDw75 (ST6GAL1) (NM_003032) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CDw75 (ST6GAL1) (NM_003032) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CDw75
Synonyms:	SIAT1; ST6Gall; ST6N
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003032
ORF Size:	1218 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208061).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003032.2
RefSeq Size:	4295 bp
RefSeq ORF:	1221 bp
Locus ID:	6480
UniProt ID:	P15907
Cytogenetics:	3q27.3
Domains:	Glyco_transf_29
Protein Families:	Secreted Protein



[View online »](#)

Protein Pathways: Metabolic pathways, N-Glycan biosynthesis

MW: 46.6 kDa

Gene Summary: This gene encodes a member of glycosyltransferase family 29. The encoded protein is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The protein, which is normally found in the Golgi but can be proteolytically processed to a soluble form, is involved in the generation of the cell-surface carbohydrate determinants and differentiation antigens HB-6, CD75, and CD76. This gene has been incorrectly referred to as CD75. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2017]