

Product datasheet for **RC208067L3V**

BPI (NM_001725) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BPI (NM_001725) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BPI
Synonyms:	BPIFD1; rBPI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001725
ORF Size:	1461 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208067).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001725.1
RefSeq Size:	1901 bp
RefSeq ORF:	1452 bp
Locus ID:	671
UniProt ID:	P17213
Cytogenetics:	20q11.23
Domains:	BPI1
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 53.9 kDa

Gene Summary: This gene encodes a lipopolysaccharide binding protein. It is associated with human neutrophil granules and has antimicrobial activity against gram-negative organisms. [provided by RefSeq, Nov 2014]