

Product datasheet for **RC208068A1V**

Human MRPL43 (NM_176792) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human MRPL43 (NM_176792) AAV Particle
Tag: Myc-DDK
Symbol: MRPL43
Synonyms: bMRP36a; L43mt; MRP-L43
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208068 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACGGCGCGGGGACTCCGAGCCGCTTCTTGCCAGCGTTCTCCACAACGGACTGGGTCGCTATGTGC
AGCAGCTGCAGCGTCTGAGCTTCAGCGTCAGCCGCGACGGCGCCTCGTCTCGCGGCGCCAGGGAGTTCGT
GGAGCGGGAGGTGATCGACTTCGCCCCGACGGAATCCAGGGGTCGTAATATATGTAACCTCGCGTCCGTGC
TGCCTGCCAGAGTAGTGCCGAATACCTTAACGGGGCTGTGCGCGAGGAGAGCATCCACTGCAAGTCGG
TCGAGGAGATCTCGACGCTGGTGCAAGCTGGCCGACCGAGTCGGGCTTGACGTGATCCGCATCCGCAA
GCCCTTCCACACCGACAACCCTAGCATCCAGGGCCAGTGGCACCCCTTCACCAACAAGCCGACCACTTC
CGCGGGCTACGCCCCGAGAGGTTAGGATCCTGCCCCAGCCAGGACACTGGCCTGAGACTGTCTGCAG
TTGACCCGAGATCCTCCTGCCCGCTGGCCCGACCCACCAGACCTCCCCACAGTGGATCCTATCTCATC
CTCATTGACCTCTGCTCCAGCCCTATGCTGTCCGAGTTTCTTGCTCCCGATTGTCCCTGCACTGACC
ACTGTGTGCTCAGCG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208068 protein sequence
Red=Cloning site Green=Tags(s)

MTARGTPSRFLASVLHNLGRYVQQLQRLSFSVSRDGAASSRGAREFVEREVIDFARRNPGVVIYVNSRPC
CVPRVVAEYLNAGVREESIHCKSVEEISTLVQKLADQSGLDVIRKPFHTDNPSIQGQWHPFTNKPTTF
RGLRPREVQDPAPAQDTGLRLSAVAPQILLPGWPDPLPTVDPISSSLTSAPAPMLSAVSCLPIVPALT
TVCSA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_176792
ORF Size:	645 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_176792.1
RefSeq Size:	2151 bp
RefSeq ORF:	648 bp
Locus ID:	84545
UniProt ID:	Q8N983
Cytogenetics:	10q24.31
MW:	23.4 kDa