

Product datasheet for **RC208082A1V**

Human DCUN1D3 (NM_173475) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DCUN1D3 (NM_173475) AAV Particle
Tag:	Myc-DDK
Symbol:	DCUN1D3
Synonyms:	44M2.4; SCCRO3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208082 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGCCAGTGTGTACCAAGTGAAGAATCCCTCATCGACCCTGGGCAGCAAAATGGAGACCGTGAGC
CCAGCAACAAGTCACATAGCAGGCGGGTGCAGGCCACCGTGAGGAGCAGGTACCACCCTGTGGCAAGCC
AGGTGGAGATATCCTCGTCAACGGGACCAAGAAGGCCGAGGCTGCCACTGAGGCCTGCCAGCTGCCAACG
TCCTCGGGAGATGCTGGGAGGGAGTCCAAGTCCAATGCCGAGGAGTCTTCCTTGCAAAGATTGGAAGAAC
TGTTCAAGCGCTACAAGGATGAGCGGGAAGATGCAATTTGGAGGAAGGCATGGAGCGCTTTTGCAATGA
CCTGTGTGTTGACCCACAGAATTTTCGAGTGCTGCTCTTGCTTGGAGTTCCAGGCTGCAACCATGTGC
AAATTCACCAGGAAGGAGTTTTTTGATGGCTGCAAAGCAATAAGTGCAGACAGCATTGACGGAATCTGTG
CACGGTTCCTAGCCTCTTAACAGAAGCCAAACAAGAGGATAAATTCAGGATCTCTACCGGTTTACATT
TCAGTTTGGCCTGGACTCTGAAGAAGGGCAGCGGTCAGTGCATCGGGAATAGCCATTGCCCTGTGGAA
CTAGTCTTTACCCAGAACATCTCCGGTATTGGACCAATGGCTAACTTCCTAACAGAGAACCCTCGG
GGATCAAGGGCATCTCCGGGACACTTGAACATGTTCTTAACCTCACTCAGGTGATTGCCCTGACCT
CAGCAACTACAGTGAAGATGAGGCCTGGCCAAGTCTTTGACACCTTTGTGGAGTGGGAAATGGAGCGA
AGGAAAAGAGAAGGGGAAGGAGAGGTGCACTCAGCTCAGGGCCTGAGGGCTTGTGTCCCGAGGAGCAGA
CT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208082 protein sequence Red=Cloning site Green=Tags(s) MGQCVTKCKNPSSSTLGSKNGDREPSNKSRRGAGHREEQVPPCGKPGGDILVNGTKKAEATEACQLPT SSGDAGRESKSNAAESSLQRLEELFRYKDEREDAILEEGMERFCNDLCVDPTEFRVLLAWKFQAATMC KFTRKEFFDGCKAISADSIDGICARFPSLLTEAKQEDKFKDLYRFTTFQGLDSEEGQSLHREIAIALWK LVFTQNNPPVLDQWLNFLTENPSGIGISRDTWNMFLNFTQVIGPDLSNYSEDEAWPSLFDTFVEWEMER RKREGEGRGALSSGPEGLCPPEQT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_173475
ORF Size:	912 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_173475.1 , NP_775746.1
RefSeq Size:	2886 bp
RefSeq ORF:	915 bp
Locus ID:	123879
UniProt ID:	Q8IWE4
Cytogenetics:	16p12.3
MW:	34.3 kDa