

## Product datasheet for **RC208101L4V**

### AP3M2 (NM\_006803) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	AP3M2 (NM_006803) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AP3M2
Synonyms:	AP47B; CLA20; P47B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006803
ORF Size:	1254 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208101).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006803.3</a>
RefSeq Size:	3624 bp
RefSeq ORF:	1257 bp
Locus ID:	10947
UniProt ID:	<a href="#">P53677</a>
Cytogenetics:	8p11.21
Domains:	Adap_comp_sub
Protein Pathways:	Lysosome



[View online »](#)

**MW:** 47 kDa

**Gene Summary:** This gene encodes a subunit of the heterotetrameric adaptor-related protein complex 3 (AP-3), which belongs to the adaptor complexes medium subunits family. The AP-3 complex plays a role in protein trafficking to lysosomes and specialized organelles. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Aug 2008]