

Product datasheet for **RC208103A1V**

Human SLAP2 (SLA2) (NM_032214) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SLAP2 (SLA2) (NM_032214) AAV Particle
Tag:	Myc-DDK
Symbol:	SLAP2
Synonyms:	C20orf156; MARS; SLAP-2; SLAP2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208103 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGAAGTCTGCCAGCAGAAGAAATCTCTGCCAAGCCCAAGCTTGAGTTCCTCTGTCCAAGGCCAGG
GACCTGTGACCATGGAAGCAGAGAGAAGCAAGGCCACAGCCGTGGCCCTGGGCAGTTTCCCGGCAGGTGG
CCCGGCCGAGCTGTCGCTGAGACTCGGGGAGCCATTGACCATCGTCTCTGAGGATGGAGACTGGTGGACG
GTGCTGTCTGAAGTCTCAGGCAGAGAGTATAACATCCCAGCGTCCACGTGGCCAAAGTCTCCCATGGGT
GGCTGTATGAGGGCTGAGCAGGGAGAAAGCAGAGGAAGTCTGTTGTTACCTGGGAACCTGGAGGGGC
CTTCCTCATCCGGGAGAGCCAGACCAGGAGAGGCTCTTACTCTGTGTCAGTCCGCCTCAGCCGCCCTGCA
TCCTGGGACCGGATCAGACACTACAGGATCCACTGCCTTGACAATGGCTGGCTGTACATCTCACC GCGCC
TCACCTTCCCCTCACTCCAGGCCCTGGTGGACCATTA CTGAGCTGGCGGATGACATCTGCTGCCTACT
CAAGGAGCCCTGTGTCCTGCAGAGGGCTGGCCCGCTCCCTGGCAAGGATATACCCCTACCTGTGACTGTG
CAGAGGACACCACTCACTGGAAGAGCTGGACAGCTCCCTCCTGTTTCTGAAGCTGCCACAGGGGAGG
AGTCTCTTCTCAGTGAGGGTCTCCGGGAGTCCCTCAGCTTCTACATCAGCCTGAATGACGAGGCTGTCTC
TTTGGATGATGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208103 protein sequence Red=Cloning site Green=Tags(s) MGSLSRRKSLPSPSLSSSVQGGPVTMEAERSKATAVALGSFPAGGPAELSLRLGEPLTIVSEDGDWWT VLSEVSGREYNIPSVHVAKVSHGWL YEGLSREKAEELLLPGNPGGAFLIRESQTRRGSYSLSVRLSRPA SWDRIRHYRIHCLDNGWLYISPRLTFPSLQALVDHYSELADDICLLKEPCVLQRAGPLPGKDIPLPVTV QRTPLNWKELDSSLLFSEAATGEESLLSEGLRESLSFYISLNDEAVSLDDA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_032214
ORF Size:	783 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032214.2</u>
RefSeq Size:	2569 bp
RefSeq ORF:	786 bp
Locus ID:	84174
UniProt ID:	<u>Q9H6Q3</u>
Cytogenetics:	20q11.23
MW:	28.6 kDa