

# Product datasheet for RC208109L3


# SENP2 (NM\_021627) Human Tagged Lenti ORF Clone

### **Product data:**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

| Product Type:                | Expression Plasmids  |
|------------------------------|--|
| Product Name:                | SENP2 (NM_021627) Human Tagged Lenti ORF Clone                 |
| Tag:                         | Myc-DDK  |
| Symbol:                      | SENP2  |
| Synonyms:                    | AXAM2; SMT3IP2   |
| Mammalian Cell<br>Selection: | Puromycin  |
| Vector:                      | pLenti-C-Myc-DDK-P2A-Puro (PS100092)                           |
| E. coli Selection:           | Chloramphenicol (34 ug/mL)                                     |
| ORF Nucleotide<br>Sequence:  | The ORF insert of this clone is exactly the same as(RC208109). |
| <b>Restriction Sites:</b>    | SgfI-Mlul  |
| Cloning Scheme:              | Cloning sites used for ORF Shuttling:                          |
|                              | GCG ATC GCC ATG NNN ACG CGT                                    |



ACCN: ORF Size:

NM\_021627 1767 bp




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

| SENP2 (NM_021627) Human Tagged Lenti ORF Clone – RC208109L3 |  |
|---|--|
| OTI Disclaimer:   | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>                                  |
| OTI Annotation:   | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| Components:   | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).   |
| Reconstitution Method:                                      | <ol> <li>Centrifuge at 5,000xg for 5min.</li> <li>Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>Close the tube and incubate for 10 minutes at room temperature.</li> <li>Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| RefSeq:   | <u>NM 021627.2</u>   |
| RefSeq Size:  | 3207 bp  |
| RefSeq ORF:   | 1770 bp  |
| Locus ID:   | 59343  |
| UniProt ID:   | <u>Q9HC62</u>  |
| Cytogenetics:   | 3q27.2   |
| Domains:  | Peptidase_C48  |
| Protein Families:   | Druggable Genome, Protease   |
| Protein Pathways:   | Wnt signaling pathway  |
| MW:   | 68.3 kDa   |
| Gene Summary:   | SUMO1 (UBL1; MIM 601912) is a small ubiquitin-like protein that can be covalently conjugated to other proteins. SENP2 is one of a group of enzymes that process newly synthesized SUMO1 into the conjugatable form and catalyze the deconjugation of SUMO1-containing species.[supplied by OMIM, Apr 2004]   |

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

# **Product images:**



Circular map for RC208109L3

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US