

## Product datasheet for **RC208134L1V**

### Mannose Phosphate Isomerase (MPI) (NM\_002435) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mannose Phosphate Isomerase (MPI) (NM_002435) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Mannose Phosphate Isomerase
Synonyms:	CDG1B; PMI; PMI1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_002435
ORF Size:	1269 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208134).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002435.1</a>
RefSeq Size:	3077 bp
RefSeq ORF:	1272 bp
Locus ID:	4351
UniProt ID:	<a href="#">P34949</a>
Cytogenetics:	15q24.1
Domains:	PMI_typel



[View online »](#)

<b>Protein Families:</b>	ES Cell Differentiation/IPS
<b>Protein Pathways:</b>	Amino sugar and nucleotide sugar metabolism, Fructose and mannose metabolism, Metabolic pathways
<b>MW:</b>	46.7 kDa
<b>Gene Summary:</b>	Phosphomannose isomerase catalyzes the interconversion of fructose-6-phosphate and mannose-6-phosphate and plays a critical role in maintaining the supply of D-mannose derivatives, which are required for most glycosylation reactions. Mutations in the MPI gene were found in patients with carbohydrate-deficient glycoprotein syndrome, type Ib. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]