

## Product datasheet for **RC208137L1V**

### JRK (NM\_003724) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	JRK (NM_003724) Human Tagged ORF Clone Lentiviral Particle
Symbol:	JRK
Synonyms:	jerky; JH8
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003724
ORF Size:	1704 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208137).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003724.2</a> , <a href="#">NP_003715.2</a>
RefSeq Size:	9116 bp
RefSeq ORF:	1707 bp
Locus ID:	8629
UniProt ID:	<a href="#">O75564</a>
Cytogenetics:	8q24.3
Domains:	DDE, CENPB, CENP-B_N
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 62.5 kDa

**Gene Summary:** This gene encodes a conserved protein that is similar to DNA-binding proteins, such as major centromere autoantigen B (CENPB). Inactivation of the related gene in mice resulted in epileptic seizures. Childhood Absence Epilepsy (CAE) has been mapped to the same chromosomal location (8q24.3) as this gene. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]