

Product datasheet for **RC208148L3V**

SERPINA12 (NM_173850) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SERPINA12 (NM_173850) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SERPINA12
Synonyms:	OL-64
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173850
ORF Size:	1242 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208148).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_173850.2
RefSeq Size:	2087 bp
RefSeq ORF:	1245 bp
Locus ID:	145264
UniProt ID:	Q8IW75
Cytogenetics:	14q32.13
Protein Families:	Druggable Genome, Secreted Protein
MW:	47.2 kDa



[View online »](#)

Gene Summary:

Adipokine that modulates insulin action by specifically inhibiting its target protease KLK7 in white adipose tissues.[UniProtKB/Swiss-Prot Function]