

## Product datasheet for **RC208154L3V**

### ETFBKMT (NM\_173802) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ETFBKMT (NM_173802) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ETFBKMT
Synonyms:	C12orf72; ETFB-KMT; METTL20
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173802
ORF Size:	786 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208154).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173802.2</a>
RefSeq Size:	2183 bp
RefSeq ORF:	789 bp
Locus ID:	254013
UniProt ID:	<a href="#">Q8IXQ9</a>
Cytogenetics:	12p11.21
Protein Families:	Druggable Genome
MW:	29.5 kDa



[View online »](#)

**Gene Summary:**

Protein-lysine methyltransferase that selectively trimethylates the flavoprotein ETFB in mitochondria (PubMed:25023281, PubMed:25416781). Thereby, may negatively regulate the function of ETFB in electron transfer from Acyl-CoA dehydrogenases to the main respiratory chain (PubMed:25416781).[UniProtKB/Swiss-Prot Function]