

## Product datasheet for **RC208165L2V**

### CCDC78 (NM\_001031737) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CCDC78 (NM_001031737) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CCDC78
Synonyms:	C16orf25; CNM4; hsCCDC78; JFP10
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001031737
ORF Size:	1314 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208165).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001031737.1</a> , <a href="#">NP_001026907.1</a>
RefSeq Size:	1611 bp
RefSeq ORF:	1317 bp
Locus ID:	124093
UniProt ID:	<a href="#">A2IDD5</a>
Cytogenetics:	16p13.3
MW:	48.5 kDa



[View online »](#)

**Gene Summary:**

The product of this gene contains two coiled-coil domains. The function of this gene is currently unknown. [provided by RefSeq, Sep 2012]