

Product datasheet for **RC208190A1V**

Human FAM78A (NM_033387) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FAM78A (NM_033387) AAV Particle
Tag:	Myc-DDK
Symbol:	FAM78A
Synonyms:	C9orf59
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208190 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGCTGTATTCAGAGCATCGGAGGCAAAGCCAGAGTCTTCCGGAAGGGATCACGGTGATTGATGTGA
AAGCCTCCATCGACCCCGTCCCCACTAGCATCGATGAGTCCTCCAGCGTGGTGCTCCGCTACCGGACACC
CCAATTCCGGGCTCGGCCAGGTGGTCATGCCGCCATCCCCAAGAAGGAGACTTGGGTAGTTGGCTGG
ATCCAGGCGTGCAGCCACATGGAGTTCTACAACCAAGTACGGCGAGCAGGGCATGTCCAGCTGGGAGCTCC
CCGACCTCCAGGAGGGCAAGATCCAAGCCATCAGCGACTCGGACGGGGTGAACCTACCCCTGGTACGGCAA
CACCACAGAGACCTGCACCATCGTGGGCCACCAAGAGGGACTCCAAGTTCATCATCAGCATGAATGAC
AACTTCTACCCAGCGTCACATGGGCCGTGCCGTGAGCGAGAGCAACGTGGCCAAGCTCACCAACATCT
ACCGGGACAGAGCTTACCACCTGGCTGGTGCCACCAACACCTCCACCAACGACATGATCATCTCTGCA
GACGCTGCACTGGCGCATGCAGCTCAGCATCGAGGTGAACCCCAACCGGCCCTGGGCCAGCGCGCCCGG
CTGCGGGAGCCCATCGCCAGGACAGCCAAAATCCTGAGCAAGAATGAGCCCATCCGCCCAGCGCCC
TGGTCAAGCCCAATGCCAACGATGCCAGGTCTCATGTGGCGGCCAAGTACGGGCAGCCGCTGGTGGT
GATCCCGCCCAAGCACCGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208190 protein sequence Red=Cloning site Green=Tags(s) MGCIQSIGGKARVFREGITVIDVKASIDPVPTSIDESSVVLRYRTPHFRASQVVMPPIPKKETWVVGW IQACSHMEFYNQYGEQGMSSWELPDLQEGKIQAISSDGVNYPWYGNTTETCTIVGPTKRDSKFIISMND NFYPSVTWAVPVSESNVAKLTNIYRDQSFTTWLVATNTSTNDMIIQLHWRMQLSIEVNPNRPLGQRAR LREPIAQDQPKILSKNEPIPPSALVKPNANDAQVLMWRPKYGQPLVVIPPKHR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_033387
ORF Size:	789 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_033387.3, NP_203745.1</u>
RefSeq Size:	3957 bp
RefSeq ORF:	852 bp
Locus ID:	286336
UniProt ID:	<u>Q5JUQ0</u>
Cytogenetics:	9q34.13
MW:	29.7 kDa