

Product datasheet for **RC208214**

TOM1L2 (NM_001033551) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TOM1L2 (NM_001033551) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TOM1L2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC208214 representing NM_001033551
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGTTCCTCTGGGAAACCGTTTCAGCACACCAGTGGGGCAGTGCCTCGAAAAGGCAACAGATGGCT
 CCCTGCAAAGTGAGGATTGGACGTTGAATATGGAGATCTGTGACATCATCAATGAGACGGAGGAAGGGCC
 AAAGGATGCCATTCGAGCCCTGAAGAAGCGGCTCAACGGGAACCGAACTACAGAGAGGTGATGCTGGCA
 TTAACAGCATGGGCTGATGCCTTTCGAAGCAGTCTGATCTACCGGCGTTGTGCACATATATGAGGAGC
 TGAAGAGGAAAGGGTTGAATTTCCCATGGCAGACTTGGACGCTCTGTCTCCCATACACACACCACAGCG
 GAGTGTCCCTGAAGTGGATCCAGCTGCGACCATGCCAGGTCCCAATCACAGCAGAGGACAAGTGTGGT
 TCCTATTCCTCGCCGCTCCTGTCCCTACTCCGACCGCAGGCCCCAGCTCTGAGTGTGACTGGCCCCA
 TCACAGCCAATTCAGAACAGATTGCCAGGCTGCGGAGTGAAGTGGACGTGTTTCGAGGAAACACAAAAGT
 CATGTCTGAGATGTTAACAGAAATGGTCCCTGGACAGGAGGATTCATCTGATCTGGAGTTGCTGCAGGAG
 CTAACAGGACCTGTCCGGCCATGCAGCAGCGCATCGTGGAGCTCATCTCCCGCGTGTCCAATGAGGAGG
 TCACCGAGGAGCTGCTGCATGTGAACGATGACCTCAACAACGTCTTCTTCGATACGAGAGGTTCGAACG
 ATACAGGTCTGGCCGATCCGTTCAAAATGCCAGTAATGGAGTACTGAATGAAGTAACCGAAGACAACCTTA
 ATAGACCTGGGGCCAGGGTCTCCAGCCGTGGTGGAGCCCAATGGTGGGGAACACAGCGCCCCCATCTTCCC
 TCTCTCCAGCTTGCAGGCTTAGACTTGGGGACAGAGAGCGTCACTGGCACCCCTCAGTTCACTCCAGCA
 ATGTAATCCCGTGACGGCTTTGACATGTTGCCAGACGAGAGGAACTCCTTGGCTGAGCAGCGCAAG
 ACGGTAACCTATGAGGATCCTCAGGCTGTCGGAGGACTTGCTTCTGCACTAGACAATCGAAAACAGAGTT
 CAGAAGGGATCCCGTTGCGCAGCCATCTGTCATGGACGACATTGAGGTGTGGCTCAGGACCGACCTGAA
 GGGTGATGATCTGGAGGAGGGTGTCAACAAGTGAAGAGTTTGATAAATTCCTTGAAGAAAGAGCCAAAGCT
 GCTGAAATGGTTCGCGACCTCCCTCGCCCCCATGGAGGCTCCTGCCCCAGCCTCAAACCTTCTGGCC
 GGAAGAAGCCAGAGCGGTGAGGATGCCCTCTTCGCCCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC208214 representing NM_001033551
 Red=Cloning site Green=Tags(s)

MEFLLGNPFSTPVGQCLEKATDGSLSQSEDWTLNMEICDIINETEEGPKDAIRALKKRLNGNRNYREVMLA
 LTAWADAFRSSPDLTGVVHIYEELKRKGVEFPMADLDALSPIHTPQRSVPEVDPAAATMPRSQSQRRTSAG
 SYSSPPPAPYSAPQAPALSVTGPITANSEQIARLRSELDVVRGNTKVMSEMLTEMVPGQEDSSDLELLQE
 LNRTCRAQQRIVELISRVSNEEVTEELLHVNDLNNVFLRYERFERYSRGRSVQNASNGVLNEVTEDNL
 IDLGPSPAVVSPMVGNTAPPSLSSQLAGLDLGTESVSGTSSLQCCNPRDGFDMFAQTRGNLSLAEQRK
 TVTYEDPQAVGGLASALDNKQSSEGIPIVAQPSVMDDIEVWLRTDLKGGDLEEGVTSEEFDKFLEERAKA
 AEMVPDLPSPMEAPAPASNPSGRKKPERSEDALFAL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_001033551

ORF Size: 1371 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033551.3](#)

RefSeq Size: 5675 bp

RefSeq ORF: 1374 bp

Locus ID: 146691


UniProt ID: [Q6ZVM7](#)

Cytogenetics: 17p11.2

MW: 49.8 kDa

Gene Summary: This gene belongs to a small gene family whose members have an N-terminal VHS domain followed by a GAT domain; domains which typically participate in vesicular trafficking. The canonical protein encoded by this gene also has a C-terminal clathrin binding motif. This protein has been shown to interact with Tollip, clathrin and ubiquitin and is thought to play a role in endosomal sorting. This gene resides in the 3.7 Mb deletion of chromosome region 17p11.2 that is associated with Smith-Magenis syndrome. Alternative splicing results in multiple transcript variants encoding distinct proteins. [provided by RefSeq, Apr 2017]

Product images:



Circular map for RC208214